



### **Product description**

Order number	55.516.177
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PP, round base, brown, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	13 mm
Height of product	100 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product brown

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 500  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 56000  
Case volume 0.0649 cbm  
Weight of product 0.0034 kg  
Weight of case 7.85 kg  
EAN of inner box 4038917108980  
EAN case 4038917031981

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

