



Product description

Order number	55.526.053
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, suitable for Abbott m1000, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Abbott m1000
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 2000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 48000
Depth of case 378 mm
Width of case 210 mm
Height of case 391 mm
Case volume 0.0309 cbm
Weight of product 0.0026 kg
Weight of case 6 kg
EAN of inner box 4038917109710
EAN case 4038917080033

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

