



### **Product description**

Order number	55.526.300
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, 500 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 500  
Type of smallest subpackaging StackPack  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 120000  
Case volume 0.0243 cbm  
Weight of product 0.0023 kg  
Weight of case 5.44 kg  
EAN of inner box 4038917367288  
EAN case 4038917367271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

