



Product description

Order number	55.528
Product description	Tube, working volume: 3 ml, (LxØ): 75 x 10 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	3 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	10 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0017 kg
Weight of case	3.89 kg
EAN of inner box	4038917082235
EAN case	4038917001328

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

