



Product description

Order number	55.598.006
Product description	Tube, working volume: 8.5 ml, (LxØ): 75 x 15.7 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8,5 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	15.7 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 56000
Case volume 0.0649 cbm
Weight of product 0.0024 kg
Weight of case 5.9 kg
EAN of inner box 4038917363136
EAN case 4038917363129

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

