



Product description

| | |
|---------------------|--|
| Order number | 57.462.015 |
| Product description | Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 250 piece(s)/StackPack |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Closure type | push cap |
| Base shape | conical base |
| Volume of work | 12 ml |
| Centrifugation max (RCF) | 2500 x g |

Size

| | |
|-------------------|--------|
| Diameter | 17 mm |
| Length of product | 110 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 42000
Depth of box 366 mm
Width of box 114 mm
Height of box 201 mm
Depth of case 388 mm
Width of case 378 mm
Height of case 250 mm
Case volume 0.0367 cbm
Weight of product 0.0063 kg
Weight of case 6.96 kg
EAN of inner box 4038917381079
EAN case 4038917381062

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

