



Product description

Order number	58.505.046
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, flat base, brown, push cap, brown, cap enclosed, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	12 ml

Size

Diameter	16.5 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	brown

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Case volume	0.0558 cbm
Weight of product	0.0064 kg
Weight of case	6.89 kg
EAN of inner box	4038917225557
EAN case	4038917225540

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

