



Product description

Order number 58.536.300
Product description Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation no
Closure type push cap
Base shape flat base
Volume of work 7 ml
Autoclavable (max. °C) 121 °C

Size

Diameter 16 mm
Height of product 50 mm
Length of product 50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 72000
Case volume 0.0492 cbm
Weight of product 0.0024 kg
Weight of case 5.24 kg
EAN of inner box 4038917367554
EAN case 4038917367547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

