



### **Product description**

Order number	58.537
Product description	Tube, working volume: 35 ml, (LxØ): 115 x 23.5 mm, material: PP, flat base, transparent, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	35 ml
Centrifugation max (RCF)	2000 x g

### **Size**

Diameter	23.5 mm
Length of product	115 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	14000
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0083 kg
Weight of case	4.78 kg
EAN of inner box	4038917082303
EAN case	4038917001403

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

