



### **Product description**

Order number	60.542.080
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, brown, screw cap, brown, cap assembled, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	16.5 mm
Length of product	57 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown

### **Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Case volume	0.0338 cbm
Weight of product	0.0046 kg
Weight of case	5.15 kg
EAN of inner box	4038917111591
EAN case	4038917079228

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

