



### Product description

Order number	60.549.001
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 66 x 11.5 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### Product characteristics

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	skirted conical base
Volume of work	3,5 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	11.5 mm
Length of product	66 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.0022 kg
Weight of case	2.81 kg
EAN of inner box	4038917112017
EAN case	4038917054928

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

