



### **Product description**

Order number	60.550.020
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, without cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	13 mm
Height of product	82 mm
Length of product	82 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 3000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 3000  
Piece(s) / pallet 48000  
Case volume 0.0813 cbm  
Weight of product 0.0029 kg  
Weight of case 10.25 kg  
EAN of inner box 4038917112284  
EAN case 4038917032537

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

