



### **Product description**

Order number	60.550.109
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	13 mm
Diameter of opening	11,55 mm
Length of product	82 mm
Length including cap	82 mm
Length excluding cap	80.9 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	40000
Depth of case	598 mm
Width of case	398 mm
Height of case	165 mm
Case volume	0.0393 cbm
Weight of product	0.0034 kg
Weight of case	3.99 kg
EAN of inner box	4038917112598
EAN case	4038917013765

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

