



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 60.551.023  |
| Product description | Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, without cap, 500 piece(s)/bag |

### **Product characteristics**

|                          |                              |
|--------------------------|------------------------------|
| Graduation               | no                           |
| Closure type             | screw cap                    |
| Base shape               | round base with skirted base |
| Volume of work           | 10 ml                        |
| Centrifugation max (RCF) | 4000 x g                     |
| Autoclavable (max. °C)   | 121 °C                       |

### **Size**

|                   |       |
|-------------------|-------|
| Diameter          | 16 mm |
| Height of product | 79 mm |
| Length of product | 79 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 2000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 32000  
Case volume 0.0813 cbm  
Weight of product 0.0024 kg  
Weight of case 5.9 kg  
EAN of inner box 4038917317801  
EAN case 4038917317795

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

