



### **Product description**

Order number	60.557.010
Product description	Screw cap tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, screw cap, without cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	4,5 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 2000  
Piece(s) / pallet 80000  
Case volume 0.0382 cbm  
Weight of product 0.0021 kg  
Weight of case 5.15 kg  
EAN of inner box 4038917113366  
EAN case 4038917056441

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

