



Product description

| | |
|---------------------|---|
| Order number | 60.9921.525 |
| Product description | Screw cap tube, working volume: 5 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, 250 piece(s)/bag |

Product characteristics

| | |
|----------------------------|-----------|
| Graduation | no |
| Base shape | flat base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 5000 x g |
| Centrifugation temperature | 21 |
| Centrifugation time | 15 Min |

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 50 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |
| Tube material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|-----------------------|
| Satisfies the requirement | sealed at 18"hg/10min |
| Latex-free | yes |
| Pyrogen-free | no |
| Safety data sheet | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 96000 |
| Depth of case | 608 mm |
| Width of case | 264 mm |
| Height of case | 324 mm |
| Case volume | 0.0519 cbm |
| Weight of case | 7 kg |
| EAN of inner box | 4038917150811 |
| EAN case | 4038917037167 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

