



### **Product description**

Order number	60.9921.807
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 97 x 16 mm, material: PP, flat base, translucent, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Base shape	flat base
Volume of work	10 ml
Centrifugation max (RCF)	5000 x g
Centrifugation temperature	21
Centrifugation time	15 Min

### **Size**

Diameter	16 mm
Length of product	97 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	translucent
Tube material	Polypropylene (PP)

### **Purity & certification**

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	608 mm
Width of case	264 mm
Height of case	324 mm
Case volume	0.0519 cbm
Weight of case	6.2 kg
EAN of inner box	4038917151009
EAN case	4038917037204

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

