



Product description

| | |
|---------------------|---|
| Order number | 60.9922.936 |
| Product description | Screw cap tube, working volume: 12 ml, (LxØ): 99 x 16 mm, material: PC, round base, transparent, screw cap, white, cap enclosed, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------------|
| Graduation | no |
| Cap | standard screw cap |
| Closure type | screw cap |
| Base shape | round base |
| Volume of work | 12 ml |
| Centrifugation max (RCF) | 8800 x g |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 99 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polycarbonate (PC) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | white |
| Tube material | Polycarbonate (PC) |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | ADR |
| Product category | no medical device IVD |
| Latex-free | yes |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 24000 |
| Depth of case | 604 mm |
| Width of case | 373 mm |
| Height of case | 264 mm |
| Case volume | 0.0595 cbm |
| Weight of product | 0.0076 kg |
| Weight of case | 8.4 kg |
| EAN of inner box | 4038917153423 |
| EAN case | 4038917027816 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

