



### **Product description**

Order number	62.476.022
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	84000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0028 kg
Weight of case	6.22 kg
EAN of inner box	4038917107402
EAN case	4038917011792

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

