



### Product description

Order number	62.493
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

### **Purity & certification**

Satisfies the requirement              sealed at 0.6 bar/10 min  
Product category                          no medical device | IVD  
Batched    no

### **Packaging**

Minimum order qty.                      2000  
Type of smallest                              bag  
subpackaging  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      2000  
Piece(s) / pallet                              72000  
Depth of case                                398 mm  
Width of case                                 398 mm  
Height of case                                265 mm  
Case volume                                 0.042 cbm  
Weight of product                         0.0031 kg  
Weight of case                                6.66 kg  
EAN of inner box                            4038917082006  
EAN case                                        4038917001113

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

