




Product Specification
Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag



Product description

Order number	62.504.037
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, with print, label/print: blue, with graduations, cap enclosed, blue, 1,000 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Label/ Print	with print
Colour of print/label	blue
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag

Size

Diameter 13 mm
Length of product 75 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap blue
Tube material Polypropylene (PP)

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 50000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 1000
Piece(s) / pallet 56000
Depth of case 568 mm
Width of case 433 mm
Height of case 120 mm
Case volume 0.0295 cbm
Weight of product 0.0038 kg
Weight of case 4.39 kg
EAN of inner box 4038917412940
EAN case 4038917412933

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

