



Product description

Order number	62.550.019
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, purple, cap enclosed, with print, label/print: purple, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	purple
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP, with print

Size

Diameter	13 mm
Length of product	82 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	purple
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Case volume	0.029 cbm
Weight of product	0.004 kg
Weight of case	4.49 kg
EAN of inner box	4038917112277
EAN case	4038917033329

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

