



**Product Specification**  
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,  
with print**

**Page 1**



**Product description**

Order number	62.551.201
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 100 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,  
with print**

**Page 2**

**Size**

Diameter	16 mm
Length of product	79 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	413 mm
Width of case	257 mm
Height of case	464 mm
Case volume	0.0492 cbm
Weight of product	0.0044 kg
Weight of case	5.04 kg
EAN of inner box	4038917112901
EAN case	4038917072113

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**