




**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**

**Page 1**



**Product description**

Order number	62.553.001
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, red, cap enclosed, with print, label/print: white/blue, with graduations, 50 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**

**Size**

Diameter	17 mm
Length of product	120 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	20000
Case volume	0.0382 cbm
Weight of product	0.0086 kg
Weight of case	4.54 kg
EAN of inner box	4038917112932
EAN case	4038917033428

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

