



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**

**Page 1**



### Product description

Order number	62.553.002
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/rack

### Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**

**Size**

Diameter 17 mm  
Length of product 120 mm

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap red

**Purity & certification**

Satisfies the requirement IATA  
Product category no medical device | IVD  
Purity standard free from DNA/DNase/RNase, free from  
pyrogens/endotoxins, non-cytotoxic, sterile  
Non-cytotoxic yes  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest rack  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 18000  
Case volume 0.0436 cbm  
Weight of product 0.0101 kg  
Weight of case 5.31 kg  
EAN of inner box 4038917112949  
EAN case 4038917033435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

