



Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**



Product description

Order number	62.554.100
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/rack

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**

Size

Diameter 17 mm
Length of product 120 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red

Purity & certification

Satisfies the requirement IATA
Product category no medical device | IVD
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation Electron irradiation
Pyrogen-free yes
Non-cytotoxic yes
Batched yes

Packaging

Minimum order qty. 25000
Type of smallest rack
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 18000
Case volume 0.0436 cbm
Weight of product 0.0074 kg
Weight of case 4.56 kg
EAN of inner box 4038917376341
EAN case 4038917376334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

