




**Product Specification**  
**Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP,  
with print**

**Page 1**



**Product description**

Order number	62.557.003
Product description	Screw cap tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, round base, brown, screw cap, brown, cap enclosed, with print, label/print: white, with graduations, 1,000 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP,  
with print**

**Page 2**

**Size**

Diameter 12 mm  
Length of product 75 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product brown  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap brown

**Purity & certification**

Satisfies the requirement IATA  
Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 1000  
Piece(s) / pallet 56000  
Case volume 0.029 cbm  
Weight of product 0.0033 kg  
Weight of case 3.76 kg  
EAN of inner box 4038917113359  
EAN case 4038917065108

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**