



### Product description

Order number	62.9924.274
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, translucent, cap, yellow, cap assembled, with printed graduations, label/print: white, with graduations, sterile, 100 piece(s)/bag

### Product characteristics

Label/ Print	with printed graduations
Colour of print/label	white
Graduation	yes
Cap	standard cap
Closure type	cap
Base shape	conical base
Volume of work	10 ml
Centrifugation max (RCF)	5000 x g
Centrifugation temperature	21
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Diameter	16 mm
Length of product	100 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	sealed at 18"hg/10min
Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	608 mm
Width of case	264 mm
Height of case	324 mm
Case volume	0.0515 cbm
Weight of case	6 kg
EAN of inner box	4038917155144
EAN case	4038917037549

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

