



Product Specification
Screw cap tube, 10 ml, (LxØ): 97 x 16 mm, PS,
with paper label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 63.9921.827 |
| Product description | Screw cap tube, working volume: 10 ml, (LxØ): 97 x 16 mm, material: PS, skirted conical base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: yellow, sterile, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|----------------------|
| Label/ Print | with paper label |
| Colour of print/label | yellow |
| Graduation | no |
| Cap | standard screw cap |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 10 ml |
| Centrifugation max (RCF) | 5000 x g |
| Centrifugation temperature | 21 |
| Centrifugation time | 15 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Screw cap tube, 10 ml, (LxØ): 97 x 16 mm, PS,
with paper label

Page 2

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 97 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |
| Tube material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|-----------------------|
| Satisfies the requirement | sealed at 18"hg/10min |
| Latex-free | yes |
| Purity standard | sterile |
| Pyrogen-free | no |
| Safety data sheet | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 48000 |
| Depth of case | 608 mm |
| Width of case | 264 mm |
| Height of case | 324 mm |
| Case volume | 0.0519 cbm |
| Weight of case | 8.4 kg |
| EAN of inner box | 4038917151108 |
| EAN case | 4038917037242 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT