




Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,
with paper label

Page 1



Product description

Order number	63.9922.248
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: white/yellow, Patient information, with graduations, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with paper label
Colour of print/label	white/yellow
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,
with paper label**

Page 2

Size

Diameter 25 mm
Length of product 90 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 300
Type of smallest blister
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 300
Piece(s) / pallet 4800
Depth of case 548 mm
Width of case 408 mm
Height of case 373 mm
Case volume 0.0834 cbm
Weight of product 0.0086 kg
Weight of case 4.3 kg
EAN of inner box 4038917238649
EAN case 4038917238632

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT