



Product description

| | |
|---------------------|---|
| Order number | 65.716.308 |
| Product description | Screw cap, violet, suitable for screw cap micro tubes, 500 piece(s)/bag |

Product characteristics

| | |
|--------------|-----------------------|
| Cap | flat cap |
| Closure type | screw cap |
| Suitable for | screw cap micro tubes |

Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | violet |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 720000 |
| EAN of inner box | 4038917925389 |
| EAN case | 4038917925396 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

