



Product description

Order number 65.716.326
Product description Screw cap, orange, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

Product characteristics

Cap flat cap
Closure type screw cap
Suitable for screw cap micro tubes

Material & colours

Product material Polypropylene (PP)
Colour of product orange
O-ring material Ethylene propylene diene rubber (EPDM)

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Type of smallest subpackaging double bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 540000
Depth of case 351 mm
Width of case 280 mm
Height of case 138 mm
Case volume 0.0136 cbm
Weight of case 2.05 kg
EAN of inner box 4038917923507
EAN case 4038917923514

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

