



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.690.378   |
| Product description | Reaction tube, working volume: 1.5 ml, material: PP, yellow, cap: yellow, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag |

### **Product characteristics**

|                          |                    |
|--------------------------|--------------------|
| Graduation               | yes                |
| Cap                      | retaining loop cap |
| Volume of work           | 1,5 ml             |
| Centrifugation max (RCF) | 20000 x g          |
| Autoclavable (max. °C)   | 121 °C             |

### **Size**

|                   |        |
|-------------------|--------|
| Sample volume     | 1.5 ml |
| Length of product | 39 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | yellow             |
| Colour of cap     | yellow             |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 160000        |
| Case volume                   | 0.0558 cbm    |
| Weight of product             | 0.001 kg      |
| Weight of case                | 5.6 kg        |
| EAN of inner box              | 4038917369534 |
| EAN case                      | 4038917369527 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

