



### **Product description**

Order number	72.699
Product description	Reaction tube, working volume: 0.5 ml, material: PP, transparent, cap: transparent, standard cap, cap attached, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard cap
Closure type	hinged lid
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	7.8 mm
Length of product	30 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	320000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0005 kg
Weight of case	2.83 kg
EAN of inner box	4038917083515
EAN case	4038917002431

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

