



Product description

| | |
|---------------------|---|
| Order number | 72.607 |
| Product description | Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, without cap, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | no |
| Knurls | yes |
| Base shape | conical base |
| Volume of work | 1,5 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 5000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 210000
Depth of case 388 mm
Width of case 378 mm
Height of case 250 mm
Case volume 0.0367 cbm
Weight of product 0.0011 kg
Weight of case 5.87 kg
EAN of inner box 4038917082839
EAN case 4038917001854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

