



### **Product description**

Order number	72.607.772
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, without cap, sterile, 500 piece(s)/double bag

### **Product characteristics**

Graduation	no
Knurls	yes
Base shape	conical base
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Satisfies the requirement IATA  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### **Packaging**

Minimum order qty. 2500  
Type of smallest subpackaging double bag  
Piece(s) / inner box 500  
Piece (s) / outer case 2500  
Piece(s) / pallet 130000  
Depth of case 548 mm  
Width of case 368 mm  
Height of case 130 mm  
Case volume 0.0262 cbm  
Weight of product 0.0011 kg  
Weight of case 3.24 kg  
EAN of inner box 4038917114288  
EAN case 4038917043809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

