



Product description

Order number	72.609.300
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 180000
Depth of case 450 mm
Width of case 450 mm
Height of case 246 mm
Case volume 0.0498 cbm
Weight of product 0.0012 kg
Weight of case 7.06 kg
EAN of inner box 4038917434287
EAN case 4038917434270

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

