



Product description

Order number	72.609.311
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, sterile, 500 piece(s)/double bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging double bag
Piece(s) / inner box 500
Piece (s) / outer case 2500
Piece(s) / pallet 140000
Case volume 0.029 cbm
Weight of product 0.0012 kg
Weight of case 3.63 kg
EAN of inner box 4038917434379
EAN case 4038917434362

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

