



### **Product description**

Order number	72.694.107
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, with printed writing space, with graduations, 500 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	retaining loop cap
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Diameter	10.8 mm
Length of product	44 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0017 kg
Weight of case	9.07 kg
EAN of inner box	4038917117654
EAN case	4038917015707

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

