



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.694.305   |
| Product description | Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap assembled, sterile, 100 piece(s)/bag |

### **Product characteristics**

|                          |                      |
|--------------------------|----------------------|
| Graduation               | no                   |
| Cap                      | standard screw cap   |
| Closure type             | screw cap            |
| Knurls                   | yes                  |
| Base shape               | skirted conical base |
| Volume of work           | 2 ml                 |
| Centrifugation max (RCF) | 20000 x g            |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 44 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |  |
|-------------------|--|
| Product material  | Polypropylene (PP)                     |
| Colour of product | transparent                            |
| Closure material  | Polypropylene (PP)                     |
| Colour of cap     | natural                                |
| O-ring material   | Ethylene propylene diene rubber (EPDM) |

### **Purity & certification**

|                           |                                |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA                           |
| Product category          | In vitro diagnostic device, CE |
| CE certified              | CE                             |
| Purity standard           | sterile                        |
| Sterilisation             | Electron irradiation           |
| Batched                   | yes                            |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 108000        |
| Depth of case                 | 351 mm        |
| Width of case                 | 280 mm        |
| Height of case                | 138 mm        |
| Case volume                   | 0.0136 cbm    |
| Weight of case                | 1.82 kg       |
| EAN of inner box              | 4038917926300 |
| EAN case                      | 4038917926317 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

