



Product description

Order number	72.730.311
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/double bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Sterilisation	Electron irradiation
Batched	yes

Packaging

Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	2500
Piece(s) / pallet	130000
Depth of case	375 mm
Width of case	545 mm
Height of case	135 mm
Case volume	0.0276 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap micro tube, 0.5 ml

Page 2

Weight of case	3.58 kg
EAN of inner box	4038917994576
EAN case	4038917994583

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT