



### **Product description**

Order number	72.785
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, transparent, without cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Knurls	no
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Satisfies the requirement IATA  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 5000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 5000  
Piece(s) / pallet 180000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 300 mm  
Case volume 0.0475 cbm  
Weight of product 0.0012 kg  
Weight of case 6.63 kg  
EAN of inner box 4038917084116  
EAN case 4038917077675

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

