



Product description

Order number	76.571.001
Descrição do artigo	Tampa de encaixe, para copo de 100 ml, 75.570.001, material: PET, transparente, 50 unid./pacote de tubo flexível

Product characteristics

Tipo de fechamento	tampa de encaixe
Indicado para	copo de 100 ml, 75.570.001

Material & colours

Material do produto	Politereftalato de etileno (PET)
Cor do produto	transparente

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tampa de encaixe, para copo de 100 ml,
75.570.001, PET, transparente

Page 2

Packaging

Quantidade mínima do pedido	1000
Tipo de embalaje interior	pacote de tubo flexível
Unidades por caixa interna	50
Unidades por caixa	1000
Unidades por palete	32000
Profundidade da embalagem externa	375 mm
Largura da embalagem externa	295 mm
Altura da embalagem externa	315 mm
Volume da embalagem externa	0.0348 cbm
Peso do produto	0.001 kg
Peso da embalagem externa	1.4 kg
EAN da caixa interna	4038917974462
EAN da embalagem externa	4038917974479

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT