



### **Product description**

Order number	65.718
Descrição do artigo	Tampa para recipiente de amostra 73.666, natural, comprimento: 8 mm, material: PEBD, 1.000 unid./pacote

### **Product characteristics**

Tampa	Tampa com saia
Tipo de fechamento	tampa de pressão
Indicado para	recipiente de amostra 73.666

### **Size**

Comprimento do produto	8 mm
------------------------	------

### **Material & colours**

Material do produto	Polietileno de baixa densidade (PEBD)
Cor do produto	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tampa para recipiente de amostra 73.666,**  
**comprimento: 8 mm, PEBD**

**Page 2**

**Purity & certification**

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

**Packaging**

Quantidade mínima do pedido	10000
Tipo de embalagem interior	pacote
Unidades por caixa interna	1000
Unidades por caixa	10000
Unidades por palete	360000
Profundidade da embalagem externa	398 mm
Largura da embalagem externa	398 mm
Altura da embalagem externa	300 mm
Volume da embalagem externa	0.0475 cbm
Peso do produto	0.0006 kg
Peso da embalagem externa	6.97 kg
EAN da caixa interna	4038917083669
EAN da embalagem externa	4038917002561

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**