



Product Specification

Tubo de rosca, 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo plano, PP, tampa incluída, 1.000 unid./pacote

Page 1



Product description

Order number	62.611.009
Descrição do artigo	Tubo de rosca, volume de trabalho: 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo plano, transparente, material: PP, com impressão, rótulo/impressão: branca, com escala, tampa incluída, natural, 1.000 unid./pacote, 1.000 unidades/caixa

Product characteristics

Rótulo/impressão	com impressão
Cor impressão/rótulo	branca
Graduação	sim
Forma do fundo	fundo falso cônico, fundo do tubo plano
Volume de trabalho	5 ml
Forma do recipiente interior	cônico

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tubo de rosca, 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo plano, PP, tampa incluída, 1.000 unid./pacote

Page 2

Size

Diâmetro	15.3 mm
Altura do produto	92 mm
Comprimento do produto	92 mm
Comprimento sem tampa	92 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	natural
Material do tubo	Polipropileno (PP)

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	1000
Unidades por caixa	1000
Unidades por palete	36000
Volume da embalagem externa	0.0492 cbm
Peso do produto	0.005 kg
Peso da embalagem externa	5.49 kg
EAN da caixa interna	4038917317894
EAN da embalagem externa	4038917317887

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT