



## Product Specification

Purga de ar para Monovette® para gases sanguíneos, transparente/laranja, comprimento: 24,5 mm

Page 1



### Product description

Order number 14.1148  
Descrição do artigo Purga de ar para Monovette® para gases sanguíneos, transparente/laranja, comprimento: 24,5 mm, acessórios para Monovette® para gases sanguíneos, 500 unid./caixa de papelão

### Product characteristics

Campo de aplicação acessórios para Monovette® para gases sanguíneos  
Exemplo de utilização Filtro para dispensa de ar residual, para gasometria Monovette®  
Indicado para Monovette® para gases sanguíneos

### Size

Comprimento do produto 24.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Purga de ar para Monovette® para gases  
sanguíneos, transparente/laranja, comprimento:  
24,5 mm**

**Page 2**

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente/laranja

### Purity & certification

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Isento de látex	sim
Em lotes	sim

### Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	500
Unidades por caixa	1000
Unidades por palete	168000
Profundidade da caixa interna	191 mm
Largura da caixa interna	129 mm
Altura da caixa interna	129 mm
Profundidade da embalagem externa	288 mm
Largura da embalagem externa	218 mm
Altura da embalagem externa	150 mm
Volume da embalagem externa	0.0094 cbm
Peso do produto	0.0009 kg
Peso da embalagem externa	1.34 kg
EAN da caixa interna	4038917085465
EAN da embalagem externa	4038917075312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**