



Product Specification
**Microvette® CB 200 NC, 200 µl, tampa preta,
tampa de pressão, fundo plano**

Page 1



Product description

Order number	18.1325.100
Descrição do artigo	Microvette® CB 200 NC, coleta de sangue capilar, preparo: Solução-tampão de tricitrato de sódio (1:4), volume nominal: 200 µl, (CxØ) com tampa: 46 x 10,8 mm, tampa preta, tampa de pressão, fundo plano, incl. pipeta de hemossedimentação, 50 unid./pacote

Product characteristics

Tipo de retirada	capilar
Tipo de preparo	Solução-tampão de tricitrato de sódio (1:4)
Código de cor	ISO
Tipo de fechamento	tampa de pressão
Campo de aplicação	coleta de sangue capilar
Forma do fundo	fundo plano

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Size

Volume nominal	200 µl
Diâmetro	10.8 mm
Comprimento do produto	40 mm
Comprimento com tampa	46 mm
Comprimento sem tampa	40 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	Low Density Polyethylene (LD-PE)
Cor da tampa	preta

Purity & certification

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

Packaging

Quantidade mínima do pedido	15000
Tipo de embalagem interior	pacote
Unidades por caixa interna	50
Unidades por caixa	500
Unidades por palete	32000
Profundidade da caixa interna	271 mm
Largura da caixa interna	126 mm
Altura da caixa interna	36 mm
Profundidade da embalagem externa	283 mm
Largura da embalagem externa	268 mm
Altura da embalagem externa	210 mm
Volume da embalagem externa	0.0159 cbm
Peso do produto	0.0036 kg
Peso da embalagem externa	2.6 kg
EAN da caixa interna	4038917127653

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® CB 200 NC, 200 µl, tampa preta,
tampa de pressão, fundo plano**

Page 3

EAN da embalagem externa 4038917021333

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT