



Product description

Order number 86.1996
Descrição do artigo Pipeta de VHS, com invólucro do êmbolo montado

Purity & certification

Categoria do produto Diagnóstico in vitro, CE
Certificação CE CE
Em lotes sim

Packaging

Quantidade mínima do pedido 1000
Tipo de embalagem interior caixa de papelão
Unidades por caixa interna 100
Unidades por caixa 1000
Unidades por palete 42000
Profundidade da caixa interna 233 mm
Largura da caixa interna 164 mm
Altura da caixa interna 76 mm
Profundidade da embalagem externa 388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipeta de VHS

Page 2

Largura da embalagem externa	338 mm
Altura da embalagem externa	250 mm
Volume da embalagem externa	0.0328 cbm
Peso do produto	0.0076 kg
Peso da embalagem externa	8.61 kg
EAN da caixa interna	4038917088138
EAN da embalagem externa	4038917006231

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT