



### **Product description**

Order number	86.1996.062
Descrição do artigo	Pipeta de VHS, com escala impressa e êmbolo montado, com impressão

### **Product characteristics**

Rótulo/impressão	com impressão
------------------	---------------

### **Purity & certification**

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

### **Packaging**

Quantidade mínima do pedido	1000
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	42000
Profundidade da caixa interna	233 mm
Largura da caixa interna	164 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipeta de VHS, com impressão**

**Page 2**

Altura da caixa interna	76 mm
Profundidade da embalagem externa	388 mm
Largura da embalagem externa	338 mm
Altura da embalagem externa	250 mm
Volume da embalagem externa	0.0328 cbm
Peso do produto	0.0076 kg
Peso da embalagem externa	8.61 kg
EAN da caixa interna	4038917135047
EAN da embalagem externa	4038917025966

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**