



Product description

Order number 90.1060.062
Descrição do artigo Rack VHS, adequado para S-Monovette®, sem traseiro fixo

Product characteristics

Versão sem traseiro fixo
Indicado para S-Monovette®

Purity & certification

Categoria do produto Diagnóstico in vitro, CE
Certificação CE CE
Em lotes sim

Packaging

Quantidade mínima do pedido 1
Unidades por caixa interna 1
Unidades por caixa 20
Unidades por palete 1560
Profundidade da embalagem externa 428 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack VHS, adequado para S-Monovette®

Page 2

Largura da embalagem externa	323 mm
Altura da embalagem externa	130 mm
Volume da embalagem externa	0.018 cbm
Peso do produto	1.167 kg
Peso da embalagem externa	23.8 kg
EAN da caixa interna	4038917124935
EAN da embalagem externa	4038917019545

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT